

Mouse Monoclonal Antibody to

## TFF3 (Trefoil Factor 3, human)

clone 15C6

CIONE	9 1200				88	- ·	+49-7641-455 670 +49-7641-455 671
Order No.:		0099-100/TFF3-15C6				Š	
Size (µg)		100				8	
Lot No.:		0099S				03/080507F	
Isotype Sp	pecies Reactivity	Applications	Mol. Weight	Ref.Cell Line	Epitope		Immunogen
lgG1 hu	uman	WB, IP, IHC (PS), ICC	9 kDa	none	N-terminus (a	a 1-12)	synthetic peptide conjugated to KLH
Background a	and Specificity:					Related Pro	ducts
ascribed to their motogenic activity in early phases of epithelial restitution. During recent years expression was found to be associated with other mucous epithelia as well. Moreover, TFF peptides appear to be involved in the enhancement of neoplastic progression. Three human genes have been identified (TFF1, TFF2 and TFF3) clusted on chromosome 21q22.3. <b>Mab TFF3-15C6</b> specifically interacts with the N-terminus of human Trefoil factor 3 (TFF3) and detects a ca. 9 kDa polypeptide in Western blots.							
Purification:	super exclu	The antibody was purified from serum-free cell culture supernatant by subsequent thiophilic adsorption and size exclusion chromatography.					
		ophilized from 1 ml 2 x PBS / 0.09 % Na-azide / PEG and acrose.					
<b>Reconstitution:</b> Reconstitute with 7			1 ml H₂O (15 min, RT).				
Stability:	Upon reconstitution, aliquote and freeze in liquid nitrogen; reconstituted antibody can be stored frozen at -80°C up to 1 year. Thaw aliquots at 37°C. Thawed aliquots may be stored at 4°C up to 3 months.						
	Avoi	Avoid repeated freeze / thaw cycles.					12
Positive Conti							
lmmunoblottii	Reco block		<mark>king buffer:</mark> Ca bation buffer, e.	isein/Tween 20 bas .g. nanoTools prodi ⊇T.		14.4	5— 4— 5—
Immunoprecip	pitation: use 1	l µg per 10⁰ SW4	80 cells			Detection of TFF-3	3
mmunocytoc	chemistry: use a	at 1 - 10 µg/ml (2	- 5 µg/ml for DA	AB staining)			were applied to SDS-PAGE and transf . The immunoblot was probed with ma
ELISA:	ND	ND					ne 2: TFF-3 overexpressed in

All products are supplied for research and investigational use only. Not for use in humans or laboratory animals.

nand TIKOERPERTECHNIK

www.nanotools.de

email: info@nanotools.de

phone: +49-7641-455 670

orders & support: