

**MEM Vitamins 100X**  
w/o L-Glutamine**CAT N°** : X0556**Theoretical pH** :  $6.5 \pm 1$  (before dilution)**Osmolality** : As reported**Colour** : Yellow**Storage conditions** : -20°C, protected from light**Shelf life** : 24 months**Sterility tests** :

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

**Endotoxin** : < 1 EU/ml**Composition** : Displayed on website ; also available on request**Recommended use** :

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light
- Manipulate the product in aseptic conditions (e.g. : under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g. : gloves, mask, hygiene cap, overall...)
- In order to preserve all product qualities, it is recommended to thaw out the flask, to aliquote, then to re-freeze the produced flasks rather than to thaw out and re-freeze the flask at each use.
- It is recommended to use the product immediately after its thaw out.

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.

**Application** :

Adding supplements of vitamins to media stimulates growth and prolongs the viability of cells in culture.

**Utilisation** :

Add this MEM Vitamins 100X w/o L-Glutamine to your basal medium (the final concentration of the MEM must be 1X) and use this medium as a classical MEM.

**Indications of deterioration** :

Medium should be clear and free of particulate and flocculent material.

Do not use if medium is cloudy or contains precipitate.

Other evidence of deterioration may include colour change or degradation of physical or performance characteristics.